


Abstract

Emory Glycomics and Molecular Interactions Core (EGMIC) is equipped with top-of-the-line LC-MS instruments. Our core provides extensive support to researchers through a broad range of mass spectrometry tools and techniques. This includes the options to learn to conduct LC-MS analysis for trainees or to have our dedicated staff run the samples. We are unique in the world by offering ion mobility mass spectrometry with electron capture dissociation. These techniques allow characterization of the shape and quaternary structure of proteins, and resolve glycan, and lipid isomers. In addition, we conduct validation of reagents especially recombinant proteins and synthetic peptides. This allows the user to confirm they received what they ordered or made the correct protein. We find 10% or more of synthetic peptides from reputable suppliers are incorrect and recombinant proteins contain surprising additions or truncations ~25% of the time. LC-MS offers precise molecular mass measurements and comprehensive analyses of proteins, peptides, glycans, and oligonucleotides. For antibodies, proteins and peptides, LC-MS provides insights into purity, folding stages, dimer configurations, protein interactions, post-translational modifications, and top-down sequencing. In glycan analysis, LC-MS enables the study of various glycoforms and glycosylation in monoclonal antibodies, utilizing workflows such as Agilent InstantPC tagging for rapid and high-throughput N-glycan composition analysis. Validation of synthetic oligonucleotides, LC-MS delivers critical data on sequence confirmation, purity, and structural modifications. Our team is dedicated to helping you leverage LC-MS to gain deeper insights into your samples, ensuring high-quality data and impactful scientific outcomes.

Application 1: Protein. Replace your Westerns, ELISA, study protein folding.


Example -1: Rapid and reproducible native protein analysis



Source: Agilent Application Note#5994-1739EN

Example -2: Native mass spectrometry to characterize intact protein and protein complex

- Native mass spectrum of the tetrameric pyruvate kinase (PK) showed two major charge envelopes. The deconvoluted spectrum revealed two multi-proteom form complex of PK tetramers: full-length PK and truncated PK tetramer.
- Native MS spectrum of Glutamate dehydrogenase (GDH). The molecular mass of the intact hexameric GDH was determined to be 334,754 Da.
- Native MS spectrum of beta-galactosidase. The molecular mass of the intact tetrameric beta-galactosidase was determined to be 465,788 Da with 1 to 10 microgram sample injection.



Source: Agilent Application Note#5994-1739EN

Application 3: Glycan, glycoform, glycopeptide in Glycosylation study

Four levels of LC/MS workflows





Figure 1. Various glycoforms/glycans quantitative analysis workflows.

Example-6: Characterization glycoforms of intact NISTmAb


- Relative quantification of mAb with 5 glycoforms
- Quantification results comparison of glycoforms from two NISTmAb samples
- Subunits analysis. Three major glycoforms were observed in heavy chain (HC)
- InstantPC labeled N-glycans from NISTmAb




Source: Agilent Application Note#5991-8796EN

Example-7: N-glycan profiling from glycoprotein using InstantPC tag labeling workflow

- High sensitivity
- Confident assignment of glycan structures
- Quick: 1 hour
- High throughput: 96 well format




Source: Agilent Application Note#5991-8071EN



Source: EGMIC Internal

Application 2: Peptide. Post-translational modifications (PTMs).


Example-3: Sensitive peptide detection



| Peptide Sequence | Mass-to-Charge Ratio |
|------------------|----------------------|
| VTSGSTTSR | 509.3 m/z |
| LASVSVSR | 428.2 m/z |
| VYVVDVAVR | 566.8 m/z |
| VYGVGVAVR | 459.8 m/z |
| LLSGAGEPK | 537.3 m/z |
| LGTLDFSK | 535.3 m/z |

Source: EGMIC Internal


Example-4: Is your peptide a dimer? Do you have a pure peptide product?



- Target peptide identified, expected mass: 4969.328
- Observed mass: Z=5, 1st isotope mass=994.8734, calculated mass: 4969.311
- This peptide is not a dimer since no presence of mass expected from a dimer (1420.808, 1105.295, 904.5142).
- This sample contained a mixture of peptides.

Example-5: Top-down amino acid sequences


- High-throughput workflow of peptide mapping provides complete amino acid sequences of mAb
- Automatic data analysis for complete sequence mapping
- Provides information on post-translational modifications (PTMs) and location



Source: Agilent Application Note#5991-7815EN


Application 4: Oligonucleotides (DNA/RNA/ASO). Impurity, truncations, additions.

Example-8: Oligonucleotides detection



Source: EGMIC Internal

Example-9: Impurities analysis



Example-10: Rapid mRNA capping analysis

- Sensitive and efficient method for process optimization and quality control of nucleic acid therapies




Figure 6. Mass spectra of ARCA Cap 0 mRNA oligonucleotides. (A - C) Deconvoluted mass spectra from Peaks 1 - 3 of Figure 5. (D - F) Deconvoluted mass spectra from Peaks 4 - 6 of Figure 5. Mass peaks marked with asterisks (16504.62 Da, 16549.69 Da, 17195.36 Da, 16564.11 Da and 17208.73 Da) were matched to their putative identities (sequences inset) with ± 20 ppm mass error (Table 5). Underlined letters in green indicate non-templated nucleotides likely due to T7 transcriptional slippage.

Source: Agilent Application Note#5994-3984EN

LC-MS Mass Spectrometry

<< Agilent LC/Q-TOF Ion Mobility 6560 IMS >>

- With ECD, CID, mass range of m/z 5-3,200
- The gold standard in direct collision cross section measurement and accuracy
- Intact proteins, native and denatured, calculate molecular size directly
- Separation of complex isobaric classes such as peptides, lipids and glycans. Used to resolve structural isomers, peptide cross-links, post-translational modifications and characterize protein shape and protein-protein interactions.

<< Agilent AdvancedBio LC/Q-TOF 6545XT >>

- High Mass Capable Q-TOF, mass range of m/z 50-30,000
- Unique ability for the top-down sequencing of proteins and synthetic peptides
- Unique fragmentation modes provides characterization of post translational modifications of peptides and structural characterization of lipids.
- Used for metabolomics, lipidomics, and qualitative flux analysis, protein-interaction mapping.

<< ThermoFisher Scientific TSQ ALTIS >>

- Triple-stage quadrupole mass spectrometer
- High sensitivity to quantify all types of molecules, such as metabolites, illicit drug, contaminants, at ultra-low levels in complex biological matrices.
- Extra selectivity for peptide quantification with Highly Selected Reaction Monitoring (H-SRM) capability.
- Used for Pharma, Biopharmaceutical, clinical research, forensic toxicology testing, food and environmental safety testing.

Molecular Interactions

<< SPR: Surface Plasmon Resonance >>



Cytiva

<< ITC: Isothermal Titration Calorimeter >>



Malvern Panalytical

<< Functional Glycan Microarray >>



Innopsys

Why us?

- Extensive support on project planning and technology consultation.
- Expert support on training for instrument operation, assay development and data interpretation.
- Meticulously maintained instruments.
- No confusion - we are here to help!
- No surprises - full transparency on costs!
- No delays - reserve and run today!