

Sample Preparation for Purine HPLC Assay

Summary

This protocol was developed to measure purines from cultured cells. Application to cultured cells was described in Lewers et al, Neuroscience, 152:761-772, 2008; Shirley et al, J Neurochem, 101, 841-853, 2007.

Sample processing

1. Sample collection
 - a. Cultured cells may be harvested by scraping or trypsinization following removal of overlying culture medium. Any excess culture medium should be removed following centrifugation at 300 × g for 5 minutes.
 - b. CSF samples should be collected and aliquoted.
 - c. Plasma samples should be collected in EDTA coated tube and centrifuged at 1000 × g for 15 minutes at 4°C to remove blood cells.
 - d. Samples may be stored indefinitely at -80°C
2. Sample deproteinization by perchloric acid
 - a. For a single 10 cm plate of cells, the volume is about 200 µL and the sample is disrupted by probe sonication.
 - b. For CSF samples, 1/10 of the sample volume of cold 1 M perchloric acid is added to sample to bring the final concentration to 0.1 M and mixed thoroughly.
 - c. For plasma samples, 1/10 of the sample volume of cold 4 M perchloric acid is added to sample to bring the final concentration to 0.4 M and mixed thoroughly.
 - d. The instructions on how to prepare perchloric acid solution:
<https://www.sigmaaldrich.com/chemistry/stockroom-reagents/learning-center/technical-library/molarity-calculator.html>
3. HPLC analysis for purines on Waters instrument
 - a. Centrifuge the homogenates at 10,000 × g for 10 minutes at 4°C.
 - b. Collect the supernatant into a new tube, being cautious to avoid collecting any of the pellet. Reserve the pellet for measurement of protein concentration.
 - c. Adjust pH of supernatant to 7.0 by adding a small volume of 2.5% 3.5 M K₂CO₃ and store on ice for 10-15 minutes to precipitate potassium perchlorate.
 - d. Centrifuge the homogenates at 10,000 × g for 10 minutes at 4°C.
 - e. Filter any remaining particulate matter from the supernatants individually by spinning in a 0.45 µm PVDF microcentrifuge filter tube at 5000 × g for 5 minutes at 4°C.
 - f. Load filtrate directly into HPLC vials for analysis on the Waters instrument.
 - g. If samples are refrozen, they must be filtered through new microcentrifuge filter tubes prior to HPLC analysis.
 - h. Part numbers of spin filters and HPLC loading vials:
0.2 µm spin filters: Catalog Number: F2517-5, Thermo Scientific. National 750 µL nonsterile micro-centrifugal PVDF membrane filters (0.2 µm pore size).
0.45 µm spin filters: Catalog Number: F2517-6, Thermo Scientific. National 750 µL nonsterile micro-centrifugal PVDF membrane filters (0.45 µm pore size).

HPLC loading vials for autosampler: Catalog Number: 186002639, Waters.
Polypropylene 12 x 32 mm Screw Neck Vial, with Cap and Preslit PTFE/Silicone
Septum, 300 μ L Volume.

4. Protein analysis
 - a. Remove all residual supernatants.
 - b. Resuspend pellet in 500 μ L of 2% SDS by gentle pipetting. Dissolving the pellet fully may require sonication for 10 seconds along with overnight incubation at 37°C.
 - c. Dilute 5-15 μ L of sample protein in additional 2% SDS to a total volume of 50 μ L for the assay.
 - d. Prepare a protein standard curve that ranges from 0 to 2 mg/ml of BSA, also in 2% SDS.
 - e. Add 200 μ L of Pierce BCA Protein Assay reagent (mix 50 parts reagent A to 1 Part of reagent B) to each 10 μ L of sample and vortex mix.
 - f. Incubate 30 minutes at 37°C.
 - g. Read on the microplate spectrophotometer at 562 nm